



## Locking Plates Systems



Bring Back To Action With 'SKY'

Locking Systems ( S. S. & Titanium)  
Locking Plates For 4mm Screws (s.s & Titanium)  
Locking Plates For 5mm Screws (s.s & Titanium)



# LOCKING SYSTEMS ( S. S. & TITANIUM)



**Locking Cortical Screw  
2.7mm (Hex 24 TPI)**  
Length : 10mm to 40 mm  
Variation 2mm  
Head : Dia. 5mm, leng. 3mm  
Drill : 2.2mm  
Cat. No. : S.S. 290.10 to 290.40  
TITA. 290T.10 to 290T.40



**Locking Cortical Screw  
3mm (Hex 24 TPI)**  
Length : 10mm to 50 mm  
Variation 2mm  
Head : Dia. 6mm, leng. 3mm  
Drill : 2.7mm  
Cat. No. : S.S. 285.10 to 285.50  
TITA. 285T.10 to 285T.50



**Locking Cortical Screw  
4mm (Hex 24 TPI)**  
Length : 10mm to 40 mm  
Variation 2mm  
Head : Dia. 6mm, leng. 3mm  
Drill : 3.2mm  
Cat. No. : S.S. 255.10 to 255.40  
TITA. 255T.10 to 255T.40



**Locking Cancellous Screw  
4mm (Hex)**  
Length : 20mm to 70 mm  
Variation 2mm  
Head : Dia. 6mm, leng. 3mm  
Drill : 3.2mm  
Cat. No. : S.S. 210.20 to 210.70  
TITA. 210T.20 to 210T.70



**Locking Cortical Screw  
5mm (Hex 24 TPI)**  
Length : 20mm to 50 mm  
Variation 2mm  
Head : Dia. 8mm, leng. 4mm  
Drill : 4mm  
Cat. No. : S.S. 256.20 to 256.50  
TITA. 256t.20 to 256T.50



**Locking Cancellous Screw  
5mm x Full th.(Hex)**  
Length : 20mm to 100 mm  
Variation 5mm  
Head : Dia. 8mm, leng. 5mm  
Drill : 4mm  
Cat. No. : S.S. 283.20 to 283.100  
TITA. 283T.20 to 283T.100



**Locking Cancellous Screw  
6.5mm x Full th.(Hex)**  
Length : 40mm to 90 mm  
Variation 5mm  
Head : Dia. 8mm, leng. 5mm  
Drill : 4.5mm  
Cat. No. : S.S. 259.40 to 259.90  
TITA. 259T.40 to 259T.90

## LOCKING PLATES FOR 4MM SCREWS (S.S & TITANIUM)



**Locking Small L.C. D. C. Plate**  
Length : 5 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3.5mm  
Cat. No. : S.S. 260.05 to 260.10  
TITA. 260T.05 to 260T.10



**Locking Reconstruction Plate**  
Length : 5 Hole to 12 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2.5mm  
Cat. No. : S.S. 261.05 to 261.12  
TITA. 261T.05 to 261T.12



**Locking Recon. Combi Hole Plate**  
Length : 4 Hole to 12 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2mm  
Cat. No. : S.S. 340.04 to 340.12  
TITA. 340T.04 to 340T.12

# LOCKING PLATES FOR 4mm SCREWS (S.S & TITANIUM)



## Locking Proximal Humerus Plate

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 14mm  
Thick. : 3mm  
Cat. No. : S.S. 264.03 to 264.10  
TITA. 264T.03 to 264T.10



## Locking Philos Plate

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 3mm  
Cat. No. : S.S. 284.03 (L-R) to 284.10 (L-R)  
TITA. 284T.03 (L-R) to 284T.10 (L-R)



## Locking Oblique 'T' Plate (Left & Right)

Length : 3 Hole to 6 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2mm  
Cat. No. : S.S. 263.03 (L-R) to 263.06 (L-R)  
TITA. 263T.03 (L-R) to 263T.06 (L-R)



## Locking Small 'T' Plate

Length : 3 Hole to 6 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2mm  
Cat. No. : S.S. 262.03 to 262.06  
TITA. 262T.03 to 262T.06



## Locking Clavicle 'S' Plate (Left & Right)

Length : 5 Hole to 12 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3mm  
Cat. No. : S.S. 288.05 (L-R) to 288.12 (L-R)  
TITA. 288T.05 (L-R) to 288T.12 (L-R)



## Locking Clavicle Medial 'S' Plate (Left & Right)

Length : 5 Hole to 12 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3mm  
Cat. No. : S.S. 339.05 (L-R) to 339.12 (L-R)  
TITA. 339T.05 (L-R) to 339T.12 (L-R)



## Locking Olecranon Plate (Left & Right)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2.5mm  
Cat. No. : S.S. 313.04 (L-R) to 313.10 (L-R)  
TITA. 313T.04 (L-R) to 313T.10 (L-R)

# LOCKING PLATES FOR 4mm SCREWS (S.S & TITANIUM)



## Locking Clavicle Hook Plate (Left & Right)

Length : 2 Hole to 6 Hole  
Variation 1 Hole  
Width : 11mm  
Thick. : 3mm  
Cat. No. : S.S. 329.02 (L-R) to 329.10 (L-R)  
TITA. 329T.02 (L-R) to 329T.10 (L-R)



## Locking Extra Articular Humerus Plate (Left & Right)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 14mm  
Thick. : 4.5mm  
Cat. No. : S.S. 332.04 (L-R) to 332.10 (L-R)  
TITA. 332T.04 (L-R) to 332T.10 (L-R)



## Locking D.H.P. Plate (Left & Right)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3mm  
Cat. No. : S.S. 278.04 (L-R) to 278.10 (L-R)  
TITA. 278T.04 (L-R) to 278T.10 (L-R)



## Locking A.L.D.H.P. Plate (Left & Right)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3mm  
Cat. No. : S.S. 277.04 (L-R) to 277.10 (L-R)  
TITA. 277T.04 (L-R) to 277T.10 (L-R)



## Locking Distal Medial Humerus Plate (L & R)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 3mm  
Cat. No. : S.S. 326.04 (L-R) to 326.10 (L-R)  
TITA. 326T.04 (L-R) to 326T.10 (L-R)



## Locking 1/3 Tubular Plate

Length : 4 Hole to 12 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2mm  
Cat. No. : S.S. 300.04 to 300.12  
TITA. 300T.04 to 300T.12



## Locking Proximal Medial Tibial "T" Plate (L & R)

Length : 4 Hole to 10 Hole  
Variation 1 Hole  
Width : 11mm  
Thick. : 4mm  
Cat. No. : S.S. 334.04 (L-R) to 334.10 (L-R)  
TITA. 334T.04 (L-R) to 334T.10 (L-R)



## Locking Low Profile Proximal Tibial Plate (L & R)

Length : 4 Hole to 12 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 4mm  
Cat. No. : S.S. 307.04 (L-R) to 307.12 (L-R)  
TITA. 307T.04 (L-R) to 307T.12 (L-R)



### Locking Postero Medial Proximal Tibia Plate

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 3.5mm  
Cat. No. : S.S. 338.03 to 338.10  
TITA. 338T.03 to 338T.10



### Locking Calcaneal "P" Plate (L & R)

Thick. : 2.5mm  
Cat. No. : S.S. 336.PL. to 336.PR  
TITA. 336T.PL. to 336T.PR



### Locking Calcaneal "Q" Plate (L & R)

Thick. : 2.5mm  
Cat. No. : S.S. 336.QL. to 336.QR  
TITA. 336T.QL. to 336T.QR



### Locking Calcaneal "R" Plate (L & R)

Thick. : 2.5mm  
Cat. No. : S.S. 336.RL. to 336.RR  
TITA. 336T.RL. to 336T.RR



### Locking Antero Lateral Distal Tibia Plate (L & R)

Length : 4 Hole to 12 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 4mm  
Cat. No. : S.S. 306.04 (L-R) to 306.12 (L-R)  
TITA. 306T.04 (L-R) to 306T.12 (L-R)



### Locking DCLC 23mm Angle Plate (L & R)

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2.5mm  
Cat. No. : S.S. 289.03 (L-R) to 289.10 (L-R)  
TITA. 289T.03 (L-R) to 289T.10 (L-R)



### Locking Distal Radius Medial Plate (L & R)

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 10mm  
Thick. : 2.5mm  
Cat. No. : S.S. 291.03 (L-R) to 291.10 (L-R)  
TITA. 291T.03 (L-R) to 291T.10 (L-R)

# LOCKING VOLAR PLATES SYSTEMS FOR 2.7MM SCREWS (S.S & TITANIUM)



## Locking Dorsal Distal Radius L Plate (L & R)

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 292.02 (L-R) to 292.05 (L-R)  
TITA. 292T.02 (L-R) to 292T.05 (L-R)



## Locking Dorsal Distal Radius T Plate

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 293.02 to 293.05  
TITA. 293T.02 to 293T.05



## Locking Dorsal Distal Radius Oblique L Plate (L & R)

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 294.02 (L-R) to 294.05 (L-R)  
TITA. 294T.02 (L-R) to 294T.05 (L-R)



## Locking Dorsal Distal Radius Right Angled L Plate (L & R)

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 295.02 (L-R) to 295.05 (L-R)  
TITA. 295T.02 (L-R) to 295T.05 (L-R)



## Locking Dorsal Distal Radius Straight Plate

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 296.02 to 296.05  
TITA. 296T.02 to 296T.05



## Locking Volar Juxta Articular Standard Plate (L & R)

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Width : 8mm  
Thick. : 2mm  
Cat. No. : S.S. 298.02 (L-R) to 298.05 (L-R)  
TITA. 298T.02 (L-R) to 298T.05 (L-R)



## Locking Volar Extra Articular Standard Plate (L & R)

Length : 2 Hole to 5 Hole  
Variation 1 Hole  
Thick. : 2mm  
Cat. No. : S.S. 297.02 (L-R) to 297.05 (L-R)  
TITA. 297T.02 (L-R) to 297T.05 (L-R)

# LOCKING PLATES FOR 5mm SCREWS (S.S & TITANIUM)



## Locking Narrow L.C.D.C. Plate

Length : 4 Hole to 14 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 4.5mm  
Cat. No. : S.S. 265.04 to 265.14  
TITA. 265T.04 to 265T.14



## Locking Broad L.C.D.C. Plate

Length : 6 Hole to 16 Hole  
Variation 1 Hole  
Width : 16mm  
Thick. : 5mm  
Cat. No. : S.S. 266.06 to 266.16  
TITA. 266T.06 to 266T.16



## Locking Reconstruction Plate

Length : 5 Hole to 12 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 2.5mm  
Cat. No. : S.S. 267.05 to 267.12  
TITA. 267T.05 to 267T.12



## Locking Recon. Combi Hole Plate

Length : 4 Hole to 12 Hole  
Variation 1 Hole  
Width : 12mm  
Thick. : 3mm  
Cat. No. : S.S. 341.04 to 341.12  
TITA. 341T.04 to 341T.12



## Locking 'T' Buttress Plate

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 16mm  
Thick. : 2.5mm  
Cat. No. : S.S. 268.04 to 268.10  
TITA. 268T.04 to 268T.10



## Locking Distal Lateral Femoral (Cobra Type) Plate (L & R)

Length : 6 Hole to 12 Hole  
Variation 1 Hole  
Width : 16mm  
Thick. : 5.5mm  
Cat. No. : S.S. 273.06 (L-R) to 273.12 (L-R)  
TITA. 273T.06 (L-R) to 273T.12 (L-R)



## Locking Cloverleaf Plate

Length : 3 Hole to 10 Hole  
Variation 1 Hole  
Width : 16mm  
Thick. : 2mm  
Cat. No. : S.S. 269.03 to 269.10



# LOCKING PLATES FOR 5mm SCREWS (S.S & TITANIUM)



## Locking Distal Lateral Tibia Plate (L & R)

Length : 4 Hole to 10 Hole  
 Width : 14mm  
 Thick. : 4mm  
 Cat. No. : S.S. 270.04 (L-R) to 270.10 (L-R)  
 TITA. 270T.04 (L-R) to 270T.10 (L-R)



## Locking Proximal Lateral Tibia Plate (L & R)

Length : 4 Hole to 12 Hole  
 Variation 1 Hole  
 Width : 14mm  
 Thick. : 4.5mm  
 Cat. No. : S.S. 271.04 (L-R) to 271.12 (L-R)  
 TITA. 271T.04 (L-R) to 271T.12 (L-R)



## Locking Condylar Buttress Plate

Length : 4 Hole to 12 Hole Left & Right  
 Variation 1 Hole  
 Width : 16mm  
 Thick. : 5mm  
 Cat. No. : S.S. 272.04 (L-R) to 272.12 (L-R)



## Locking Proximal Femoral Plate

Length : 4 Hole to 14 Hole  
 Variation 1 Hole  
 Width : 17mm  
 Thick. : 6mm  
 Cat. No. : S.S. 274.04 to 274.14  
 TITA. 274T.04 to 274T.14



## Locking 'L' Buttress Plate (L & R)

Length : 3 Hole to 12 Hole  
 Variation 1 Hole  
 Width : 16mm  
 Thick. : 2.5mm  
 Cat. No. : S.S. 276.03 (L-R) to 276.12 (L-R)  
 TITA. 276T.03 (L-R) to 276T.12 (L-R)



## Locking Distal Lateral Femoral (New Type) (L & R)

Length : 4 Hole to 12 Hole  
 Variation 1 Hole  
 Width : 18mm  
 Thick. : 6mm  
 Cat. No. : S.S. 333.04 (L-R) to 333.12 (L-R)  
 TITA. 333T.04 (L-R) to 333T.12 (L-R)



## Locking D.H.S. Barrel Plate

Angle : 1250, 1300, 1350, 140  
 Length : 4 Hole to 14 Hole  
 Variation 1 Hole  
 Barrel Length : 32mm Dai. : 12.5mm  
 Width : 19mm  
 Thick. : 6mm  
 Cat. No. : S.S. 156.125-04 to 156.140-14  
 TITA. 156T.125-04 to 156T.140-14



Email: [skysurgicals50@yahoo.in](mailto:skysurgicals50@yahoo.in) / [office@skysurgicals.com](mailto:office@skysurgicals.com)  
Web.: [www.skysurgicals.com](http://www.skysurgicals.com)