



## Nailing Systems



Bring Back To Action With 'SKY'

Interlocking Nail & Bolt For Tibia & Femur  
Instruments For Tibia & Femur Nail  
Interlocking Nail & Bolt For Supracondylar  
Instruments For Supracondylar Nail  
Interlocking Nail & Bolt For T.f./p.f. & Long P.f.  
Instruments For T.f./p.f. & Long P.f. Nail  
Interlocking Nail & Bolt For Humerus  
Instruments For Humerus Nail  
Ender Nailing Systems  
Instruments For Enders Nail  
Elastic Nailing Systems  
Instruments For Elastic Nail





**Inter Locking Nail For Tibia**

Dia : 8, 9, 10 & 11mm  
 Length : 28cm to 38cm  
 Variation 1cm  
 Cat. No. : S.S. 161.828 to 161.1138  
 Tita. 161T.828 to 161T.1138



**Inter Locking Nail For Femur**

Dia : 9, 10 & 11mm  
 Length : 34cm to 44cm  
 Variation 1cm  
 Cat. No. : S.S. 162.934 to 161.1144  
 Tita. 162T.934 to 161T.1144



**Inter Locking Nail For Tibia (High Band)**

Dia : 8, 9 & 10mm  
 Length : 28cm to 38cm  
 Variation 1cm  
 Cat. No. : S.S. 161H.828 to 161H.1038  
 Tita. 161HT.828 to 161HT.1038



**I. L. Bolt 3.9mm (Hex)**

Length : 24cm to 50cm  
 Variation 2mm  
 Head : Dia. : 8mm, Leng. : 4mm  
 Cat. No. : S.S. 163.24 to 163.50  
 Tita. 163T.24 to 163T.50



**I. L. Bolt 4.9mm (Hex)**

Length : 28mm to 60mm  
 Variation 2mm  
 Head : Dia. : 8mm, Leng. : 4mm  
 Cat. No. : S.S. 164.28 to 164.60  
 Tita. 164T.28 to 164T.60



**End Cap For Tibia & Femur Nail**

Cat. No. : S.S.168.01S  
 Tita. 168.01T

# INSTRUMENTS FOR TIBIA & FEMUR NAIL



**Introducer with Th. bolt for Tibia & Femur**  
Cat. No. : 411.01



**Extractor Rod**  
Cat. No. : 411.03



**Jig for Proximal Locking Tibia & Femur**  
Cat. No. : 411.02



**Jig For Distal locking Tibia**  
Cat. No. : 411.17



**Round Hammer**  
Cat. No. : 411.04



**Jig for Distal locking Tibia High Band**  
Cat. No. : 403.03



**Protection Sleeve**  
Cat. No. : 411.05



**Jig for Distal Locking Femur**  
Cat. No. : 411.18



**Drill Sleeve for 3.2mm Drill**  
Cat. No. : 411.06



**Pointer**  
Cat. No. : 411.08



**Drill Sleeve for 4mm Drill**  
Cat. No. : 411.07



**Tissue protector**  
Cat. No. : 409.07



**Spanner ordinary (2 pcs.)**  
Cat. No. : 411.12

## INSTRUMENTS FOR TIBIA & FEMUR NAIL



**Screw Driver 4.5mm (Long)**  
Cat. No. : 411.13



**Depth Gauge Long**  
Cat. No. : 411.09



**Guide Wire 30" Long** 2mm & 3mm  
Cat. No. : 411.14-15



**Bone Awl Diamond Pointed**  
Cat. No. : 409.06



**Box for Instruments**  
Cat. No. : 411.16



**Drill Bits 10" Long**  
3.2 mm & 4 mm  
Cat. No. : 401.2732-40

## INTERLOCKING NAIL & BOLT FOR SUPRACONDYLAR



**Supracondylar Nail**  
Dia. : 9, 10 & 11  
Length : 15, 20, 25 & 30 cm  
Variation : 5cm  
Cat. No. : 167.0915 to 167.1130



**I.L. Bolt 4.9mm (Hex)**  
Length : 28mm to 60mm  
Variation : 2mm  
Head Dia. : 8mm, length 4mm  
Cat. No. : 164.28 to 164.60

## INSTRUMENTS FOR SUPRACONDYLAR NAIL



**Distal & Proximal Locking Jig Left**  
Cat. No. : 413.01L



**Distal & Proximal Locking Jig Right**  
Cat. No. : 413.01R



# INSTRUMENTS FOR SUPRACONDYLAR NAIL



**Introducer with Thr. Bolt**  
Cat. No. : 413.02



**Round Hammer**  
Cat. No. : 412.04



**Extractor Rod**  
Cat. No. : 412.03



**Screw Driver 4.5mm Long**  
Cat. No. : 411.13



**Pointer**  
Cat. No. : 411.08



**Spanner Ordinary (2 Pcs.)**  
Cat. No. : 412.09



**Drill Sleeve For 4mm Drill**  
Cat. No. : 411.07



**Drill Bits 4mm x 10" Long**  
Cat. No. : 401.2740



**Protection Sleeve**  
Cat. No. : 411.05



**Guide Wire 3mm**  
Cat. No. : 411.15



**Awl Diamond Pointed**  
Cat. No. : 409.06



**Box For Instruments**  
Cat. No. : 413.10



**Depth Gauge Long**  
Cat. No. : 411.09

# INTERLOCKING NAIL & BOLT FOR T.F./P.F. & LONG P.F.



## Trochanteric Femoral Nail

Dia. : 9,10,11 & 12mm  
 Length : 180mm  
 Angle : 130° & 135°  
 Cat. No. : S.S. 275.130-09 to 275.135-12  
 Tita. 275T.130-09 to 275T.135-12



## Cannulated Cervical 8 mm Screw

Length : 70mm to 110mm  
 Variation : 5 mm  
 Head Dia. : 10mm len. 5mm  
 Cat. No. : S.S. 232.70 to 232.110  
 Tita. 232T.70 to 232T.110



## Proximal Femoral Nail

Dia. : 9,10,11 & 12mm  
 Length : 250mm  
 Angle : 130° & 135°  
 Cat. No. : S.S. 230.130-09 to 230.135-12  
 Tita. 230T.130-09 to 230T.135-12



## Cannulated Stabilization 6.4 mm Screw

Length : 60mm to 100mm  
 Variation : 5 mm  
 Head Dia. : 8mm len. 5mm  
 Cat. No. : S.S. 233.60 to 233.100  
 Tita. 233T.60 to 233T.100



## Long Proximal Femoral Nail

Left & Right  
 Dia. : 9,10,11 & 12mm  
 Length : 360mm to 420mm Variation 20mm  
 Angle : 130° & 135°  
 Cat. No. : 231L-R. 130-0936 TO 231L-R. 135-1242



## Distal locking Bolt for 4.9mm (Hex) (I.L. Bolt)

Length : 28mm to 60mm  
 Variation : 2 mm  
 Head Dia. : 8mm len. 4mm  
 Cat. No. : S.S. 164.28 to 164.60  
 Tita. 164T.28 to 164T.60



## End Cap For T.F. & P.F. NAIL

Cat. No. : S.S.168.02S  
 TITA 168.02T

# INSTRUMENTS FOR T.F./P.F. & LONG P.F. NAIL



**P.F.N. Jig With Distal 130°**  
Cat. No. : 424.01130



**P.F.N. Jig With Distal 135°**  
Cat. No. : 424.01135



**Protection Sleeve For Cervicals 8mm Screw for 10 mm**  
Cat. No. : 424.03



**Introducer With Thr. Bolt**  
Cat. No. : 424.02



**Protection Sleeve For Cervicals 8mm Screw for 8 mm**  
Cat. No. : 424.04



**Extractor Rod**  
Cat. No. : 424.08



**Protection Sleeve For 6.4mm Screw for 8 mm**  
Cat. No. : 424.05



**Round Hammer With Handle**  
Cat. No. : 411.04



**Protection Sleeve For Distal Locking Screw 8 mm**  
Cat. No. : 411.05



**Step Reamer For 8mm Screw**  
Cat. No. : 424.09



**Drill Sleeve For Distal Locking Screw 4mm**  
Cat. No. : 411.07



**Step Reamer For 6.4mm Screw**  
Cat. No. : 424.10



**Guide Wire Sleeve For Cervical Screw 2.5mm**  
Cat. No. : 424.06



**Guide Wire Sleeve For Stabilization Screw 2.mm**  
Cat. No. : 424.07



# INSTRUMENTS FOR T.F./P.F. & LONG P.F. NAIL



**Cannulated Tap 6.4mm & 8mm**  
Cat. No. : 424.11-12



**Hammering Device Set of 3 Pc**  
Cat. No. : 424.01



**Cannulated Screw Driver For 8mm & 6.4mm Screw**  
Cat. No. : 424.13



**Entry Reamer 15mm**  
Cat. No. : 424.14



**'T' Handle Quick Coupling For Reamer**  
Cat. No. : 424.19



**Awl Diamond Pointed Cannulated**  
Cat. No. : 424.15



**Screw Driver Long 4.5mm**  
Cat. No. : 411.13



**Drill Bits 4mm x 10" Long for 4.9mm Screw**  
Cat. No. : 401.2740



**Depth Gauge Long**  
Cat. No. : 424.23



**Guide Wire 16" Long** (Size : 2mm, 2.5mm)  
Cat. No. : 424.16-17



**Solapuri Sleeve**  
Cat. No. : 424.21



**Guide Wire 30" Long** (Size : 3mm)  
Cat. No. : 411.15



**T.F.N. Jig With Distal Locking 130°**  
Size : 424.24130



**T.F.N. Jig With Distal Locking 135°**  
Size : 424.24135

## INSTRUMENTS FOR T.F./P.F. & LONG P.F. NAIL



**Spanner Ordinary 2nos.**

Cat. No. : 411.12



**'T' Handle Spanner For Introducer**

Cat. No. : 424.18



**Box For Instruments**

Cat. No. : 424.22

## INTERLOCKING NAIL & BOLT FOR HUMERUS



**Inter locking Nail for Humerus**

Dia : 6.5, 7 & 7.5mm

Length : 20 cm to 30 cm

Variation : 2cm

Cat. No. : S.S. 165.6520 to 165.7530

Tita. 165T.6520 to 165T.7530



**I.L. Bolt for 3.9mm (Hex)**

Length : 24mm to 50mm

Variation : 2mm

Head Dia : 6mm, leng 4mm

Cat. No. : S.S. 166.24 to 166.50

Tita. 166T.24 to 166T.50



**I.L. Bolt for 2.9mm (Hex)**

Length : 24mm to 50mm

Variation : 2mm

Head Dia : 6mm, leng 3mm

Cat. No. : S.S. 287.24 to 287.50

Tita. 287T.24 to 287T.50

# INSTRUMENTS FOR HUMERUS NAIL



**Jig for Proximal locking**  
Cat. No. : 412.01



**Pointer**  
Cat. No. : 411.08



**Introducer with Threaded Bolt**  
Cat. No. : 412.02



**Drill Sleeve 3.2mm**  
Cat. No. : 411.06



**Extractor Rod**  
Cat. No. : 412.03



**Drill Sleeve 2.5mm**  
Cat. No. : 411.19



**Round Hammer**  
Cat. No. : 412.04



**Awl for Humerus 8mm**  
Cat. No. : 412.06



**Protection Sleeve**  
Cat. No. : 411.05



**Depth Gauge Long**  
Cat. No. : 411.09



**Jig for Distal Locking Humerus**  
Cat. No. : 412.12



**Drill Bits Extra Long 8" 3.2mm**  
Cat. No. : 401.2632

## INSTRUMENTS FOR HUMERUS NAIL



**Guide Wire 24" Long 2mm**  
Cat. No. : 412.10



**Spanner Ordinary (2 pcs.)**  
Cat. No. : 412.09



**Box for Instruments**  
Cat. No. : 412.11



**Screw driver 3.5mm Long**  
Cat. No. : 412.08

## ENDER NAILING SYSTEMS



### Ender Nail Tibia, Femur & Humerus

Dia : 2.5mm, 3.5mm & 4.5mm  
Length : 20cm to 44cm  
Variation : 1cm  
Cat. No. : 160.2520 to 160.4544



## INSTRUMENTS FOR ENDERS NAIL



**Ender Nail Extractor**  
Cat. No. : 410.0125-35



**Ender Nail Bender Pair**  
Cat. No. : 410.05



**Ender Nail Introducer Small Tips**  
Cat. No. : 410.03



**Ender Nail Awl Diamond - Pointed**  
Cat. No. : 410.06



**Ender Nail Introducer Large Tips**  
Cat. No. : 410.04

# ELASTIC NAILING SYSTEMS



## Elastic Nail Titanium 45cm Long

Dia : 2mm to 4.5mm  
Variation : 0.5 mm  
Cat. No. : 252.20 to 252.45

## INSTRUMENTS FOR ELASTIC NAIL



### Extractor Road With Threaded Bolt & Tommy Bar

Cat. No. : 422.01



### Round Hammer With Handel

Cat. No. : 422.04



### Reduction Instrument

Cat. No. : 422.02



### Bone Awl

Cat. No. : 422.06



### "T" Handle With Chuck

Cat. No. : 422.07



### Introducer Punch (Size : 2.5mm, 3.5mm, 4.5mm)

Cat. No. : 422.05-25-35-45



### Extractor Pliers Type

Cat. No. : 422.03



### Final Punch

Cat. No. : 422.08





Email: [skysurgicals50@yahoo.in](mailto:skysurgicals50@yahoo.in) / [office@skysurgicals.com](mailto:office@skysurgicals.com)  
Web.: [www.skysurgicals.com](http://www.skysurgicals.com)