



Spinel Implants



Bring Back To Action With 'SKY'

**Pedicular System For Spine (s. S. & Titanium)**  
**Instrument For Pedicular System**



# PEDICULAR SYSTEM FOR SPINE (S. S. & TITANIUM)



## Mono Axial Pedicle Screw (S.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 236.530 S to 236.650 S  
**Tita.** : 244.530 S to 244.650 S



## Mono Axial Pedicle Screw (D.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 236.530 D to 236.650 D  
**Tita.** : 244.530 D to 244.650 D



## Poly Axial Pedicle Screw (S.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 237.530 S to 237.650 S  
**Tita.** : 245.530 S to 245.650 S



## Poly Axial Pedicle Screw (D.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 237.530 D to 237.650 D  
**Tita.** : 245.530 D to 245.650 D



## Reduction Axial Pedicle Screw (S.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 235.530 S to 235.650 S  
**Tita.** : 243.530 S to 243.650 S



## Reduction Axial Pedicle Screw (D.L.)

Length : 30mm to 50mm  
Variation 5mm  
Dia. : 5mm & 6mm  
Cat. No. : **S.S.** : 235.530 D to 236.650 D  
**Tita.** : 243.530 D to 243.650 D

# PEDICULAR SYSTEM FOR SPINE (S. S. & TITANIUM)



## Connecting Rod Plain

Length : 80mm to 260mm  
Variation 20mm  
Dia. : 5mm  
Cat. No. : **S.S.** : 238.80 to 238.260  
**Tita.** : 246.80 to 246.260



## Connecting Rod Serrated

Length : 80mm to 260mm  
Variation 20mm  
Dia. : 5mm  
Cat. No. : **S.S.** : 239.80 to 239.260  
**Tita.** : 247.80 to 247.260



## Hook Wide Blade

Size : Facet, Supra laminar, Infra laminar  
Cat. No. : **S.S.** : 240.1 to 240.3  
**Tita.** : 248.1 to 248.3

## Transverse Connector

Cat. No. : **S.S.** : 241  
**Tita.** : 249



## Spinal Cage (Titanium)

Length : 30mm to 50mm  
Dia : 14mm & 16mm  
Variation 10mm  
Cat. No. : 251.1430 to 251.1650



## Connecting Rod

Length : 50mm to 80mm  
Variation 10mm  
Dia. : 3mm  
Cat. No. : **S.S.** : 242.50 to 242.80  
**Tita.** : 250.50 to 250.80

# INSTRUMENT FOR PEDICULAR SYSTEM



**Starter Awl**

Cat. No. : 425.01



**Hook Holding Forceps**

Cat. No. : 425.08



**Pedicle Centralizer 4Pc.**

Cat. No. : 425.02



**Introducer For Mono Axial Screw (S.L.)**

Cat. No. : 425.09



**Probe Round and Straight**

Cat. No. : 425.03



**Introducer For Poly Axial and Reduction Screw (S.L.)**

Cat. No. : 425.10



**Probe Flat and Curved**

Cat. No. : 425.04



**Self Holding Screw Driver For Inny & Outy**

Cat. No. : 425.11



**Bone Tap 5mm & 6mm**

Cat. No. : 425.05-06



**Screw Driver For Final Fixation of Inny Fixation**

Cat. No. : 425.12



**Pedicle Sound**

Cat. No. : 425.07



**Reduction Screw Breaker Pair**

Cat. No. : 425.13

# INSTRUMENT FOR PEDICULAR SYSTEM



**Rod Bender Roller Type**  
Cat. No. : 425.14



**Sagittal Approximator**  
Cat. No. : 425.22



**Rod Holding Forceps For 3mm & 5mm ROD**  
Cat. No. : 425.15-16



**Sagittal Rod Approximator**  
Cat. No. : 425.23



**Vice Grip**  
Cat. No. : 425.17



**Screw Distractor**  
Cat. No. : 425.24



**Rod Bender Pair**  
Cat. No. : 425.18



**Screw Compressor**  
Cat. No. : 425.25



**Rod Stabilizer**  
Cat. No. : 425.19



**Hook Inserter For Single Lock**  
Cat. No. : 425.26



**Rod Pusher Straight & Curved**  
Cat. No. : 425.20-21



**Hook Pusher**  
Cat. No. : 425.27

# INSTRUMENT FOR PEDICULAR SYSTEM



**Facet Finder**

Cat. No. : 425.28



**Poly & Reduction Screw Introducer For Double Lock**

Cat. No. : 425.34



**Supra Lamina Locator**

Cat. No. : 425.29



**Introducer For Mono Axial Screw (D.L.)**

Cat. No. : 425.35



**Infra Lamina Locator**

Cat. No. : 425.30



**Final Fixer For Outy (D.L.)**

Cat. No. : 425.36



**Staple Washer Inserter**

Cat. No. : 425.31



**Final Tightener For Double Lock**

Cat. No. : 425.37



**Rod Cutter**

Cat. No. : 425.32



**Hook Inserter For Double Lock**

Cat. No. : 425.38



**Introducer For Mono Axial Screw (D.L.)**

Cat. No. : 425.33



**Instruments Box**

Cat. No. : 425.39



**Pedicular Screw Box with Tray**

Cat. No. : 425.40



Email: [skysurgicals50@yahoo.in](mailto:skysurgicals50@yahoo.in) / [office@skysurgicals.com](mailto:office@skysurgicals.com)  
Web.: [www.skysurgicals.com](http://www.skysurgicals.com)